

SOFRELINER Precautions

The following precautions will enhance the longevity, durability and safety of your denture lined with SOFRELINER material:

1. SOFRELINER material should not be scrubbed with a denture brush or toothbrush. The best way to clean SOFRELINER is with cotton or a soft piece of gauze. While running cold water on the SOFRELINER surface, wipe the surface with the cotton or gauze dipped in soap. Use a gentle wiping motion to remove all debris and food deposits from the SOFRELINER material. The tooth surface and hard acrylic denture base may be cleaned with a regular denture brush by scrubbing with soap, denture toothpaste or regular denture cleaners.
2. Soak the denture once a day, preferably at bedtime, in an oxygenating denture cleanser (i.e., Polident, Efferdent, etc.). The denture should not be cleaned with oxygenating cleanser within the patient's mouth. Color fading may occur with long term use, depending upon the patient's denture maintenance.
3. Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER, but it can be used continuously.
4. The denture should be stored in water when not being worn.
5. The lined denture should be inspected and ultrasonically cleaned in our office every three months, or when plaque deposit is noted.
6. If plaque is routinely present on the denture, the denture should be soaked at least every other day, in plaque-control mouthwash.
7. In the event of any sign of sensitivity to silicone, use of the denture should be discontinued, and the patient should contact our office immediately.